



Product description

Order number	80.1303
Product description	Swabs, in the tube, no medium, (LxØ): 175 x 12 mm, stem material: aluminium, stem length: 133 mm, swab material: viscose, sterile, 100 piece(s)/case

Product characteristics

Closure type	push cap
Design	in the tube
Type of medium	no medium

Size

Diameter	12 mm
Length of product	175 mm
Length including cap	175 mm
Length excluding cap	155 mm
Length including cap and plunger	164 mm
Length of stem	133 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	orange
Swab material	viscose
Stem material	aluminium

Purity & certification

Product category	Medical device (MD)
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	no
Batched	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	case
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	32000
Depth of box	193 mm
Width of box	115 mm
Height of box	138 mm
Depth of case	590 mm
Width of case	300 mm
Height of case	210 mm
Case volume	0.0372 cbm
Weight of product	0.0042 kg
Weight of case	5.42 kg
EAN of inner box	4038917086110
EAN case	4038917234191

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

