



## Product description

Order number	80.1362
Product description	Swabs, with tube included, Amies gel with charcoal, (LxØ): 175 x 12 mm, stem material: PS, stem length: 131 mm, swab material: viscose, sterile, 1 piece(s)/blister

## Product characteristics

Closure type	push cap
Design	with tube included
Type of medium	Amies gel with charcoal

## Size

Diameter	12 mm
Length of product	175 mm
Length including cap	175 mm
Length excluding cap	155 mm
Length including cap and plunger	166 mm
Length of stem	131 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material	Polypropylene (PP)
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	black
Swab material	viscose
Stem material	Polystyrene (PS)

### **Purity & certification**

Product category	Medical device
CE certified	CE supplier
Purity standard	sterile
Sterilisation	Gamma irradiation
Safety data sheet	yes
Batched	yes

### **Packaging**

Minimum order qty.	500
Type of smallest subpackaging	blister
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	11500
Depth of case	520 mm
Width of case	450 mm
Height of case	220 mm
Case volume	0.0515 cbm
Weight of product	0.0113 kg
Weight of case	6.62 kg
EAN of inner box	4038917086332
EAN case	4038917004893

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

