



Product description

Order number	80.625
Product description	Swabs, in the tube, no medium, (LxØ): 108 x 16 mm, stem material: PS, stem length: 83 mm, swab material: viscose, sterile, 500 piece(s)/bag

Product characteristics

Type of cap	screw cap
Design	in the tube
Type of medium	no medium

Size

Diameter	16 mm
Length of product	108 mm
Length including cap	108 mm
Length excluding cap	101 mm
Length including cap and plunger	97 mm
Length of stem	83 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Cap material	Polyethylene (PE)
Colour of cap	transparent
Swab material	viscose
Stem material	Polystyrene (PS)

Purity & certification

Product category	In-vitro diagnostic, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	500
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	500
Piece(s) / pallet	22000
Depth of case	568 mm
Width of case	348 mm
Height of case	155 mm
Case volume	0,0306 cbm
Weight of product	0,0064 kg
Weight of case	3,73 kg
EAN of inner box	4038917083003
EAN case	4038917001984

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

