



**Product description**

Order number	92.180.610
Product description	Plate rotor, for 23 tubes Ø 8?12 mm, for SARMIX® M 2000, 1 piece(s)/case

**Product characteristics**

Suitable for	SARMIX® M 2000
Storage capacity	23

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Plate rotor, for 23 tubes Ø 8?12 mm, for SARMIX®**  
**M 2000**

Case volume	0.0122 cbm
Weight of product	0.1809 kg
Weight of case	2.01 kg
EAN of inner box	4038917101981
EAN case	4038917075749

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

