



Product Specification

Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M2000

Page 1



Product description

Order number	92.180.611
Product description	Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M2000

Product characteristics

Suitable for	SARMIX® M2000
Storage capacity	2

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	10
Case volume	0,0122 cbm
Weight of product	0,2009 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

Plate rotor, for 2 tubes up to 35 mm Ø and 6 tubes, up to 20 mm Ø and 6 tubes, up to 12.5 mm Ø, for SARMIX® M2000

Page 2

Weight of case	2,21 kg
EAN of inner box	4038917101998
EAN case	4038917075756

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-01-20

This document was prepared by EDP support and is valid without signature.



SARSTEDT