



Product description

| | |
|---------------------|---|
| Order number | 92.180.612 |
| Product description | Block rotor, for 40 tubes up to 11.5 mm diameter, for SARMIX® M 2000, 1 piece(s)/case |

Product characteristics

| | |
|------------------|----------------|
| Suitable for | SARMIX® M 2000 |
| Storage capacity | 40 |

Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Batched | no |

Packaging

| | |
|-------------------------------|------|
| Minimum order qty. | 1 |
| Type of smallest subpackaging | case |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 10 |
| Piece(s) / pallet | 10 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Block rotor, for 40 tubes up to 11.5 mm diameter,
for SARMIX® M 2000

| | |
|-------------------|---------------|
| Depth of case | 350 mm |
| Width of case | 240 mm |
| Height of case | 150 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.1509 kg |
| Weight of case | 1.71 kg |
| EAN of inner box | 4038917102001 |
| EAN case | 4038917075763 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

