



### **Product description**

Order number	92.180.613
Product description	Block rotor, for 24 tubes up to 15 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

### **Product characteristics**

Suitable for	SARMIX® M 2000
Storage capacity	24

### **Purity & certification**

Product category	no medical device   IVD
Batched	no

### **Packaging**

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**Block rotor, for 24 tubes up to 15 mm diameter,**  
**for SARMIX® M 2000**

Case volume	0.0122 cbm
Weight of product	0.1509 kg
Weight of case	1.71 kg
EAN of inner box	4038917102018
EAN case	4038917075770

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

