



Product description

Order number	92.180.614
Product description	Block rotor, for 7 tubes up to 28 mm diameter, for SARMIX® M 2000, 1 piece(s)/case

Product characteristics

Suitable for	SARMIX® M 2000
Storage capacity	7

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Block rotor, for 7 tubes up to 28 mm diameter, for
SARMIX® M 2000

Case volume	0.0122 cbm
Weight of product	0.2009 kg
Weight of case	2.21 kg
EAN of inner box	4038917102025
EAN case	4038917075787

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

