



Product description

Order number	92.180.617
Product description	Block rotor, for 24 tubes up to 17 mm Ø (15 ml tubes), for SARMIX® M 2000, 1 piece(s)/case

Product characteristics

Suitable for	SARMIX® M 2000
Storage capacity	24

Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	10
Piece(s) / pallet	10

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Block rotor, for 24 tubes up to 17 mm Ø (15 ml tubes), for SARMIX® M 2000

Page 2

Case volume	0.0122 cbm
Weight of product	0.2299 kg
Weight of case	2.31 kg
EAN of inner box	4038917245906
EAN case	4038917245890

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-01-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT