




Product Specification
Rotor, for 36 collection tubes, for micro tube Ø 10.8 mm

Page 1



Product description

Order number 92.181.610
Product description Rotor, for 36 collection tubes, for micro tube Ø 10.8 mm

Product characteristics

Suitable for micro tube Ø 10.8 mm
Storage capacity 36

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece (s) / outer case 1
Piece(s) / pallet 40
Case volume 0.0056 cbm
Weight of product 1.3595 kg
Weight of case 1.36 kg

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Rotor, for 36 collection tubes, for micro tube Ø
10.8 mm

Page 2

EAN of inner box 4038917156233
EAN case 4038917156226

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT