



Product description

Order number 92.181.612
Product description Rotor, for 6 collection tubes, for tubes Ø 28 mm

Product characteristics

Suitable for tubes Ø 28 mm
Storage capacity 6

Purity & certification

Product category no medical device | IVD
Batched yes

Packaging

Minimum order qty. 1
Piece(s) / inner box 1
Piece(s) / case 1
Piece(s) / pallet 40
Case volume 0.0053 cbm
Weight of product 1.3595 kg
Weight of case 1.36 kg
EAN of inner box 4038917156271

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
Rotor, for 6 collection tubes, for tubes Ø 28 mm

Page 2

EAN case

4038917156264

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2022-07-07

This document was prepared by EDP support and is valid without signature.



SARSTEDT