



### **Product description**

|                     |                |
|---------------------|----------------|
| Order number        | 90.181.601     |
| Product description | Rotation mixer |

### **Purity & certification**

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Batched          | yes                     |

### **Packaging**

|                        |               |
|------------------------|---------------|
| Minimum order qty.     | 1             |
| Piece(s) / inner box   | 1             |
| Piece (s) / outer case | 1             |
| Piece(s) / pallet      | 1             |
| Weight of product      | 0.0005 kg     |
| EAN of inner box       | 4038917156219 |
| EAN case               | 4038917156202 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

