



### Product description

Order number	90.181.601
Product description	Rotation mixer

### Purity & certification

Batched	yes
---------	-----

### Packaging

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Weight of product	0.0005 kg
EAN of inner box	4038917156219
EAN case	4038917156202

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

