



Product description

Order number	95.983
Product description	IsoFreeze® MCT Rack, material: PP, multi-coloured, format: 6 x 4, suitable for reaction tubes Ø 10.8 mm (1.5 ml & 2 ml), 1 piece(s)/case

Product characteristics

Label/ Print	with print
Suitable for	reaction tubes Ø 10.8 mm (1.5 ml & 2 ml)
Format	6 x 4

Material & colours

Product material	Polypropylene (PP)
Colour of product	multi-coloured
Colour of cap	transparent
Colour of frame	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product category	no medical device IVD
Batched	no

Packaging

Minimum order qty.	1
Type of smallest subpackaging	case
Piece(s) / inner box	1
Piece (s) / outer case	5
Piece(s) / pallet	240
Depth of case	197 mm
Width of case	93 mm
Height of case	584 mm
Case volume	0.0107 cbm
Weight of product	0.349 kg
Weight of case	2.5 kg
EAN of inner box	4038917394895
EAN case	4038917394888

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

