



Product Specification

IsoFreeze® PCR Rack, PP, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates

Page 1



Product description

Order number	95.984
Product description	IsoFreeze® PCR Rack, material: PP, multi-coloured, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates, Colour change from purple to pink at 7 °C, 2 piece(s)/case

Product characteristics

Suitable for	0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates
Format	12 x 8

Material & colours

Product material	Polypropylene (PP)
Colour of product	multi-coloured
Colour of cap	transparent
Colour of frame	white

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

IsoFreeze® PCR Rack, PP, format: 12 x 8, suitable for 0.1 ml and 0.2 ml PCR tubes, PCR strips and PCR plates

Page 2

Purity & certification

Product category no medical device | IVD
Batched no

Packaging

Minimum order qty. 2
Type of smallest subpackaging case
Piece(s) / inner box 2
Piece (s) / outer case 10
Piece(s) / pallet 1280
Depth of case 197 mm
Width of case 175 mm
Height of case 584 mm
Case volume 0.0201 cbm
Weight of product 0.2144 kg
Weight of case 2.95 kg
EAN of inner box 4038917373425
EAN case 4038917373418

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-05-21

This document was prepared by EDP support and is valid without signature.



SARSTEDT