



### **Product description**

Order number 95.988.002  
Product description Base station, material: PC, blue, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

### **Product characteristics**

Suitable for work tray 95.987.002 & 0.2 ml PCR tubes

### **Size**

Width of product 85.5 mm  
Height of product 30.4 mm  
Length of product 128 mm

### **Material & colours**

Product material Polycarbonate (PC)  
Colour of product blue

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Purity & certification**

Product category                      no medical device | IVD  
Batched                                      no

### **Packaging**

Minimum order qty.                      5  
Type of smallest                              case  
subpackaging  
Piece(s) / inner box                      5  
Piece (s) / outer case                      75  
Piece(s) / pallet                              3600  
Depth of case                              598 mm  
Width of case                                398 mm  
Height of case                               145 mm  
Case volume                                0.0345 cbm  
Weight of product                        0.0302 kg  
Weight of case                              2.83 kg  
EAN of inner box                        4038917124171  
EAN case                                    4038917019132

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

