



## Product description

Order number	95.988
Product description	Base station, material: PC, transparent, (LxWxH): 128 x 85.5 x 30.4 mm, suitable for work tray 95.987.002 & 0.2 ml PCR tubes, 5 piece(s)/case

## Product characteristics

Suitable for	work tray 95.987.002 & 0.2 ml PCR tubes
--------------	---

## Size

Width of product	85.5 mm
Height of product	30.4 mm
Length of product	128 mm

## Material & colours

Product material	Polycarbonate (PC)
Colour of product	transparent

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**

**Base station, PC, suitable for work tray 95.987.002  
& 0.2 ml PCR tubes**

**Purity & certification**

Product category                      no medical device | IVD  
Batched                                    no

**Packaging**

Minimum order qty.                      5  
Type of smallest                            case  
subpackaging  
Piece(s) / inner box                      5  
Piece (s) / outer case                    75  
Piece(s) / pallet                          3600  
Depth of case                              598 mm  
Width of case                               398 mm  
Height of case                              145 mm  
Case volume                                0.0345 cbm  
Weight of product                        0.0302 kg  
Weight of case                             2.83 kg  
EAN of inner box                          4038917084864  
EAN case                                    4038917003636

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

