



Product description

Order number	65.989
Product description	PCR strip of lids, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), Biosphere® plus, material: PP, highly transparent, flat cap, suitable for qPCR, 12 piece(s)/bag

Product characteristics

Cap	flat cap
Closure type	push cap
Suitable for	PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992)
Suitable for qPCR	yes
Number of wells	8
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR strip of lids, Biosphere® plus, PP, highly transparent, flat cap

Page 2

Size

Diameter	4.47 mm
Width of product	7.3 mm
Height of product	3.25 mm
Length of product	77.2 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	highly transparent

Purity & certification

Product category	In-vitro diagnostic (IVD)
CE certified	CE - manufacturer's self-declaration
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	240
Type of smallest subpackaging	bag
Piece(s) / inner box	12
Piece (s) / outer case	240
Piece(s) / pallet	77760
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0003 kg
Weight of case	0.26 kg
EAN of inner box	4038917281485
EAN case	4038917281478

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

