



Product description

Order number 65.989.002
Product description PCR strip of lids, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), PCR Performance Tested, material: PP, highly transparent, flat cap, suitable for qPCR, 120 piece(s)/minigrip bag

Product characteristics

Cap flat cap
Closure type push cap
Suitable for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992)

Suitable for qPCR yes
Number of wells 8
Autoclavable (max. °C) 121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Diameter	4.47 mm
Width of product	7.3 mm
Height of product	3.25 mm
Length of product	77.2 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	highly transparent

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	155520
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0004 kg
Weight of case	0.38 kg
EAN of inner box	4038917259910
EAN case	4038917259903

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

