



## **Product description**

Order number	72.735.002
Product description	PCR single tube, 0.5 ml, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 500 piece(s)/minigrip bag

## **Product characteristics**

Graduation	yes
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	1
Volume of work	25 µl - 0,5 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR single tube, 0.5 ml, PCR Performance**  
**Tested, transparent, PP, flat cap**

**Page 2**

**Size**

Sample volume	0.5 ml
Diameter	7.9 mm
Width of product	10 mm
Height of product	29.5 mm
Length of product	30.9 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	192000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0004 kg
Weight of case	1.1 kg
EAN of inner box	4038917119986
EAN case	4038917016988

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**