



### Product description

Order number	72.735.002
Product description	PCR single tube, 0.5 ml, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, optimised for qPCR, 500 piece(s)/bag

### Product characteristics

Graduation	yes
Cap	flat cap
Type of cap	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High-Profile
Number of wells	1
Volume of work	25 µl - 0,5 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**PCR single tube, 0.5 ml, PCR Performance**  
**Tested, transparent, PP, flat cap**

**Page 2**

**Size**

Sample volume	0,5 ml
Diameter	7,9 mm
Width of product	10 mm
Height of product	29,5 mm
Length of product	30,9 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	IVD product
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of inner packaging	bag
Piece(s) / inner box	500
Piece(s) / case	2000
Piece(s) / pallet	192000
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0,0135 cbm
Weight of product	0,0004 kg
Weight of case	1,1 kg
EAN of inner box	4038917119986
EAN case	4038917016988

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-04

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**