



**Product Specification**  
**PCR single tube, 0.5 ml, Biosphere® plus,**  
**transparent, PP, flat cap**

**Page 1**



**Product description**

Order number	72.735.100
Product description	PCR single tube, 0.5 ml, Biosphere® plus, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 100 piece(s)/bag

**Product characteristics**

Graduation	yes
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	1
Volume of work	25 µl - 0,5 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR single tube, 0.5 ml, Biosphere® plus,**  
**transparent, PP, flat cap**

**Page 2**

**Size**

Sample volume	0.5 ml
Diameter	7.9 mm
Width of product	10 mm
Height of product	29.5 mm
Length of product	30.9 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	500
Piece (s) / outer case	1000
Piece(s) / pallet	144000
Depth of box	201 mm
Width of box	151 mm
Height of box	126 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0004 kg
Weight of case	0.75 kg
EAN of inner box	4038917120050
EAN case	4038917056465

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

