



Product Specification
PCR single tube, 0.2 ml, PCR Performance
Tested, transparent, PP, flat cap

Page 1



Product description

Order number	72.737.002
Product description	PCR single tube, 0.2 ml, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 500 piece(s)/minigrip bag

Product characteristics

Graduation	no
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	1
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
PCR single tube, 0.2 ml, PCR Performance
Tested, transparent, PP, flat cap

Page 2

Size

Sample volume	0.2 ml
Diameter	7.5 mm
Width of product	7.5 mm
Height of product	20.65 mm
Length of product	21.7 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	500
Piece (s) / outer case	2000
Piece(s) / pallet	648000
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.0002 kg
Weight of case	0.53 kg
EAN of inner box	4038917120104
EAN case	4038917017060

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2024-03-28

This document was prepared by EDP support and is valid without signature.



SARSTEDT