



Product description

| | |
|---------------------|--|
| Order number | 72.737.002 |
| Product description | PCR single tube, 0.2 ml, PCR Performance Tested, transparent, material: PP, with anti-contamination shield, flat cap, optimised for qPCR, 500 piece(s)/bag |

Product characteristics

| | |
|----------------------------|----------------|
| Graduation | no |
| Cap | flat cap |
| Type of cap | push cap |
| Suitable for qPCR | yes |
| Base shape | conical |
| High or low profile | High-Profile |
| Number of wells | 1 |
| Volume of work | 10 µl - 0,2 ml |
| Centrifugation max (RCF) | 8000 x g |
| Centrifugation temperature | 22 °C |
| Centrifugation time | 5 Min |
| Autoclavable (max. °C) | 121 °C |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
PCR single tube, 0.2 ml, PCR Performance
Tested, transparent, PP, flat cap

Page 2

Size

| | |
|-------------------|----------|
| Sample volume | 0,2 ml |
| Diameter | 7,5 mm |
| Width of product | 7,5 mm |
| Height of product | 20,65 mm |
| Length of product | 21,7 mm |

Material & colours

| | |
|-------------------|--------------------|
| Product material | Polypropylene (PP) |
| Colour of product | transparent |

Purity & certification

| | |
|------------------|------------------------|
| Product category | IVD product |
| CE certified | CE |
| Purity standard | PCR Performance Tested |
| Batched | yes |

Packaging

| | |
|-------------------------|---------------|
| Minimum order qty. | 2000 |
| Type of inner packaging | bag |
| Piece(s) / inner box | 500 |
| Piece(s) / case | 2000 |
| Piece(s) / pallet | 648000 |
| Depth of case | 268 mm |
| Width of case | 191 mm |
| Height of case | 100 mm |
| Case volume | 0,0051 cbm |
| Weight of product | 0,0002 kg |
| Weight of case | 0,53 kg |
| EAN of inner box | 4038917120104 |
| EAN case | 4038917017060 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT