



**Product Specification**  
**PCR single tube, 0.2 ml, Biosphere® plus,**  
**transparent, PP, flat cap**

**Page 1**



### Product description

Order number	72.737
Product description	PCR single tube, 0.2 ml, Biosphere® plus, transparent, material: PP, with anti-contamination shield, flat cap, suitable for qPCR, 250 piece(s)/bag

### Product characteristics

Graduation	no
Cap	flat cap
Closure type	snap-on lid
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	1
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**PCR single tube, 0.2 ml, Biosphere® plus,**  
**transparent, PP, flat cap**

**Page 2**

**Size**

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	7.5 mm
Width of product	7.5 mm
Height of product	20.65 mm
Length of product	21.7 mm

**Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	2000
Piece(s) / pallet	288000
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0.01 cbm
Weight of product	0.0002 kg
Weight of case	0.72 kg
EAN of inner box	4038917083836

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**PCR single tube, 0.2 ml, Biosphere® plus,**  
**transparent, PP, flat cap**

**Page 3**

EAN case

4038917028370

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-27

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**