



Product Specification
PCR single tube, 0.2 ml, Biosphere® plus,
transparent, PP, flat cap

Page 1



Product description

Order number	72.737
Product description	PCR single tube, 0.2 ml, Biosphere® plus, transparent, material: PP, with anti-contamination shield, flat cap, optimised for qPCR, 250 piece(s)/bag

Product characteristics

Graduation	no
Cap	flat cap
Type of cap	Snap-on lid
Suitable for qPCR	yes
Base shape	conical
High or low profile	High-Profile
Number of wells	1
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
PCR single tube, 0.2 ml, Biosphere® plus,
transparent, PP, flat cap

Page 2

Size

Sample volume	0,2 ml
Max. volume	200 µl
Diameter	7,5 mm
Width of product	7,5 mm
Height of product	20,65 mm
Length of product	21,7 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	IVD product
CE certified	CE
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Batched	yes

Packaging

Minimum order qty.	2000
Type of inner packaging	bag
Piece(s) / inner box	1000
Piece(s) / case	2000
Piece(s) / pallet	288000
Depth of box	189 mm
Width of box	151 mm
Height of box	97 mm
Depth of case	323 mm
Width of case	213 mm
Height of case	145 mm
Case volume	0,01 cbm
Weight of product	0,0002 kg
Weight of case	0,72 kg
EAN of inner box	4038917083836
EAN case	4038917028370

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-09

This document was prepared by EDP support and is valid without signature.



SARSTEDT