



### Product description

Order number	72.982.002
Product description	PCR strip of 8, up to 100 µl, PCR Performance Tested, transparent, material: PP, flat lid, optimised for qPCR, 125 piece(s)/bag

### Product characteristics

Graduation	no
Cap	flat lid
Type of cap	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	Low-Profile
Number of wells	8
Volume of work	5 µl - 0,1 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22 °C
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Sample volume	0,1 ml
Max. volume	100 µl
Width of product	16,05 mm
Length of product	76,2 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

### Packaging

Minimum order qty.	125
Type of inner packaging	bag
Piece(s) / inner box	125
Piece(s) / case	1250
Piece(s) / pallet	56250
Case volume	0,0299 cbm
Weight of product	0,0016 kg
Weight of case	3,18 kg
EAN of inner box	4038917395465
EAN case	4038917395458

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

