



Product description

Order number	72.985.992
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, red/green/blue/violet, material: PP, without cap, suitable for qPCR, per colour, 120 piece(s)/minigrip bag

Product characteristics

Graduation	no
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	8
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	5.4 mm
Width of product	8.1 mm
Height of product	20.7 mm
Length of product	70.9 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	colour mix

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	155520
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0.0051 cbm
Weight of product	0.001 kg
Weight of case	0.64 kg
EAN of inner box	4038917124119
EAN case	4038917067683

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

