



Product description

Order number	72.985
Product description	PCR plate without skirt, 12 strips of 8 tubes (72.985.002) preloaded in PCR work tray (95.987.002), 96 well, transparent, High Profile, 200 µl, Biosphere® plus, material: PP, 1 piece(s)/blister

Product characteristics

Closure type	push cap
Design	12 strips of 8 tubes (72.985.002) preloaded in PCR work tray (95.987.002)
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	without skirt
High or low profile	High Profile
Number of wells	96
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification

PCR plate without skirt, 96 well, transparent, High Profile, 200 µl, Biosphere® plus, PP

Page 2

Centrifugation time 5 Min
Autoclavable (max. °C) 121 °C

Size

Sample volume 0.2 ml
Max. volume 200 µl
Diameter 5.4 mm
Width of product 79 mm
Height of product 20.7 mm
Length of product 115 mm

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard Biosphere® plus
Sterilisation Gas sterilisation
Pyrogen-free yes
Batched yes

Packaging

Minimum order qty. 20
Type of smallest subpackaging blister
Piece(s) / inner box 20
Piece (s) / outer case 20
Piece(s) / pallet 3360
Depth of case 288 mm
Width of case 218 mm
Height of case 150 mm
Case volume 0.0094 cbm
Weight of product 0.0228 kg
Weight of case 0.76 kg
EAN of inner box 4038917084840

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**PCR plate without skirt, 96 well, transparent, High
Profile, 200 µl, Biosphere® plus, PP**

Page 3

EAN case

4038917003612

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-04-26

This document was prepared by EDP support and is valid without signature.



SARSTEDT