



### **Product description**

Order number	72.991.102
Product description	PCR strip of 8, High-Profile, PCR Performance Tested, transparent, material: PP, flat cap

### **Product characteristics**

Graduation	no
Cap	flat cap
Type of cap	Snap-on lid
Base shape	conical
High or low profile	High-Profile
Number of wells	8
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	20 °C
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**PCR strip of 8, PCR Performance Tested,**  
**transparent, PP, flat cap**

**Page 2**

**Size**

Sample volume 0,2 ml

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard PCR Performance Tested  
Batched yes

**Packaging**

Minimum order qty. 1200  
Piece(s) / inner box 120  
Piece(s) / case 1200  
Piece(s) / pallet 40320  
Case volume 0,0354 cbm  
Weight of product 0,0018 kg  
Weight of case 2,77 kg  
EAN of inner box 4038917231657  
EAN case 4038917231640

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:  
2021-12-09

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**