



**Product description**

Order number	72.991.992
Product description	PCR strip of 8, up to 200 µl, High Profile, PCR Performance Tested, red/green/blue/violet, material: PP, with anti-contamination shield, flat cap, cap attached, suitable for qPCR, 120 piece(s)/minigrip bag

**Product characteristics**

Graduation	no
Cap	flat cap
Closure type	push cap
Suitable for qPCR	yes
Base shape	conical
High or low profile	High Profile
Number of wells	8
Volume of work	10 µl - 0,2 ml
Centrifugation max (RCF)	8000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	7.5 mm
Width of product	10.2 mm
Height of product	20.65 mm
Length of product	72.4 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	colour mix

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	480
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	480
Piece (s) / outer case	480
Piece(s) / pallet	34560
Depth of case	348 mm
Width of case	278 mm
Height of case	140 mm
Case volume	0.0135 cbm
Weight of product	0.0015 kg
Weight of case	1.01 kg
EAN of inner box	4038917124317
EAN case	4038917077613

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

