



Product description

Order number	72.994.002
Product description	PCR strip of 8, High Profile, PCR Performance Tested, transparent, material: PP, round cap, cap attached, 120 piece(s)/bag

Product characteristics

Graduation	no
Cap	round cap
Closure type	push cap
Base shape	conical
High or low profile	High Profile
Number of wells	8

Size

Sample volume	0.2 ml
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
PCR strip of 8, PCR Performance Tested,
transparent, PP, round cap

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 240
Type of smallest subpackaging bag
Piece(s) / inner box 240
Piece (s) / outer case 1200
Piece(s) / pallet 1200
Case volume 0.0299 cbm
Weight of product 0.0031 kg
Weight of case 4.3 kg
EAN of inner box 4038917124324
EAN case 4038917080736

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2025-04-29

This document was prepared by EDP support and is valid without signature.



SARSTEDT