



### **Product description**

Order number	72.994.002
Product description	PCR strip of 8, High Profile, PCR Performance Tested, transparent, material: PP, round cap, cap attached, 120 piece(s)/bag

### **Product characteristics**

Graduation	no
Cap	round cap
Closure type	push cap
Base shape	conical
High or low profile	High Profile
Number of wells	8

### **Size**

Sample volume	0.2 ml
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**PCR strip of 8, PCR Performance Tested,**  
**transparent, PP, round cap**

**Page 2**

**Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

**Purity & certification**

Product category no medical device | IVD  
Purity standard PCR Performance Tested  
Batched yes

**Packaging**

Minimum order qty. 240  
Type of smallest bag  
subpackaging  
Piece(s) / inner box 240  
Piece (s) / outer case 1200  
Piece(s) / pallet 1200  
Case volume 0.0299 cbm  
Weight of product 0.0031 kg  
Weight of case 4.3 kg  
EAN of inner box 4038917124324  
EAN case 4038917080736

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-21

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**