



Product description

Order number	72.994.002
Product description	PCR strip of 8, PCR Performance Tested, transparent, material: PP, round cap

Product characteristics

Graduation	no
Cap	round cap
Type of cap	push cap
Base shape	conical
High or low profile	High-Profile
Number of wells	8

Size

Sample volume	0,2 ml
---------------	--------

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Product Specification
PCR strip of 8, PCR Performance Tested,
transparent, PP, round cap

Page 2

Material & colours

Product material Polypropylene (PP)
Colour of product transparent

Purity & certification

Product category no medical device | IVD
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 240
Piece(s) / inner box 240
Piece(s) / case 1200
Case volume 0,0299 cbm
Weight of product 0,0031 kg
Weight of case 4,3 kg
EAN of inner box 4038917124324
EAN case 4038917080736

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

Date of issue:
2021-12-05

This document was prepared by EDP support and is valid without signature.



SARSTEDT