



Product description

| | |
|---------------------|---|
| Order number | 72.994.992 |
| Product description | PCR strip of 8, High Profile, PCR Performance Tested, red/green/blue/violet/yellow, material: PP, round cap, cap attached, 120 piece(s)/bag |

Product characteristics

| | |
|---------------------|--------------|
| Graduation | no |
| Cap | round cap |
| Closure type | push cap |
| Base shape | conical |
| High or low profile | High Profile |
| Number of wells | 8 |

Size

| | |
|---------------|--------|
| Sample volume | 0.2 ml |
|---------------|--------|

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material Polypropylene (PP)
Colour of product colour mix

Purity & certification

Product category no medical device | IVD
Purity standard PCR Performance Tested
Batched yes

Packaging

Minimum order qty. 1200
Type of smallest subpackaging bag
Piece(s) / inner box 120
Piece (s) / outer case 1200
Piece(s) / pallet 1200
Case volume 0.0299 cbm
Weight of product 0.0031 kg
Weight of case 4.3 kg
EAN of inner box 4038917124331
EAN case 4038917080743

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

