



### Product description

Order number 72.994.992  
Product description PCR strip of 8, PCR Performance Tested, colour mix, material: PP, round cap

### Product characteristics

Graduation no  
Cap round cap  
Type of cap push cap  
Base shape conical  
High or low profile High-Profile  
Number of wells 8

### Size

Sample volume 0,2 ml

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product colour mix

### **Purity & certification**

Product category no medical device | IVD  
Purity standard PCR Performance Tested  
Batched yes

### **Packaging**

Minimum order qty. 1200  
Piece(s) / inner box 120  
Piece(s) / case 1200  
Case volume 0,0299 cbm  
Weight of product 0,0031 kg  
Weight of case 4,3 kg  
EAN of inner box 4038917124331  
EAN case 4038917080743

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

