



Product description

Order number	65.1998.400
Product description	PCR strip of lids, for PCR plates (REF 72.1979.102/.132, 72.1982.202), PCR Performance Tested, material: PP, transparent, flat cap, suitable for qPCR, 12 piece(s)/minigrip bag

Product characteristics

Cap	flat cap
Closure type	push cap
Suitable for	PCR plates (REF 72.1979.102/.132, 72.1982.202)
Suitable for qPCR	yes
Number of wells	8
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Width of product	7.4 mm
Length of product	80.5 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	1200
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	300
Piece (s) / outer case	1200
Piece(s) / pallet	151200
Depth of case	355 mm
Width of case	255 mm
Height of case	123 mm
Case volume	0.0111 cbm
Weight of product	0.0004 kg
Weight of case	0.69 kg
EAN of inner box	4038917336369
EAN case	4038917336352

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

