



Product description

Order number	72.1977.202
Product description	PCR plate without skirt, 96 well, transparent, Low Profile, 100 µl, PCR Performance Tested, material: PP, 20 piece(s)/minigrip bag

Product characteristics

Closure type	film
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	without skirt
High or low profile	Low Profile
Number of wells	96
Volume of work	5 µl - 0,1 ml
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	0.1 ml
Max. volume	100 µl
Width of product	77.7 mm
Length of product	117.9 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	100
Type of inner packaging	minigrip bag
Piece(s) / inner box	20
Piece(s) / case	100
Piece(s) / pallet	5000
Case volume	0.0096 cbm
Weight of product	0.0128 kg
Weight of case	1.5 kg
EAN of inner box	4038917405102
EAN case	4038917405096

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

