




## Product Specification

**PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 72.1977.232  |
| Product description | PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, material: PP, optimised for qPCR, 20 piece(s)/minigrip bag |

### Product characteristics

|                                 |               |
|---------------------------------|---------------|
| Closure type                    | film          |
| Suitable for qPCR               | yes           |
| Base shape                      | conical       |
| Construction of the plate edges | without skirt |
| High or low profile             | Low Profile   |
| Number of wells                 | 96            |
| Volume of work                  | 5 µl - 0,1 ml |
| Centrifugation max (RCF)        | 3000 x g      |
| Centrifugation temperature      | 22            |
| Centrifugation time             | 5 Min         |
| Autoclavable (max. °C)          | 121 °C        |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**PCR plate without skirt, 96 well, white, Low Profile, 100 µl, PCR Performance Tested, PP**

**Page 2**

### Size

|                   |          |
|-------------------|----------|
| Sample volume     | 0.1 ml   |
| Max. volume       | 100 µl   |
| Width of product  | 77.7 mm  |
| Length of product | 117.7 mm |

### Material & colours

|                   |                    |
|-------------------|--------------------|
| Product material  | Polypropylene (PP) |
| Colour of product | white              |

### Purity & certification

|                  |                         |
|------------------|-------------------------|
| Product category | no medical device   IVD |
| Purity standard  | PCR Performance Tested  |
| Batched          | yes                     |

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 100           |
| Type of smallest subpackaging | minigrip bag  |
| Piece(s) / inner box          | 20            |
| Piece (s) / outer case        | 100           |
| Piece(s) / pallet             | 5000          |
| Depth of case                 | 137 mm        |
| Width of case                 | 137 mm        |
| Height of case                | 493 mm        |
| Case volume                   | 0.0093 cbm    |
| Weight of product             | 0.0106 kg     |
| Weight of case                | 1.5 kg        |
| EAN of inner box              | 4038917405133 |
| EAN case                      | 4038917405126 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

