



### Product description

Order number	72.1979.201
Product description	PCR plate half skirt, 96 well, transparent, High Profile, 200 µl, Biosphere® plus, material: PP, 1 piece(s)/blister

### Product characteristics

Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	half skirt
High or low profile	High Profile
Number of wells	96
Volume of work	10 µl - 200 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### Size

Sample volume	0.2 ml
Max. volume	200 µl
Diameter	4.95 mm
Width of product	86 mm
Height of product	19.33 mm
Length of product	126.1 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

### Purity & certification

Product category	no medical device   IVD
Purity standard	Biosphere® plus
Sterilisation	Gas sterilisation
Batched	yes

### Packaging

Minimum order qty.	20
Type of inner packaging	blister
Piece(s) / inner box	10
Piece(s) / case	20
Piece(s) / pallet	1120
Depth of case	282 mm
Width of case	280 mm
Height of case	247 mm
Case volume	0.0195 cbm
Weight of product	0.0169 kg
Weight of case	0.97 kg
EAN of inner box	4038917572958
EAN case	4038917572941

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

