



### **Product description**

Order number	72.1981
Product description	PCR plate half skirt, 96 well, transparent, Low Profile, 100 µl, PCR Performance Tested, material: PP, 10 piece(s)/minigrip bag

### **Product characteristics**

Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	half skirt
High or low profile	Low Profile
Number of wells	96
Volume of work	5 µl - 100 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Size

Sample volume	0.1 ml
Max. volume	100 µl
Diameter	5.48 mm
Width of product	85.4 mm
Height of product	17.37 mm
Length of product	127.8 mm

## Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

## Purity & certification

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

## Packaging

Minimum order qty.	100
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	50
Piece (s) / outer case	100
Piece(s) / pallet	9800
Depth of box	270 mm
Width of box	230 mm
Height of box	96 mm
Depth of case	282 mm
Width of case	210 mm
Height of case	247 mm
Case volume	0.0146 cbm
Weight of product	0.014 kg
Weight of case	1.75 kg
EAN of inner box	4038917573474
EAN case	4038917573467

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

