



**Product description**

Order number	72.1984.202
Product description	PCR plate full skirt, 384 well, transparent, Low Profile, 40 µl, PCR Performance Tested, material: PP, 50 piece(s)/minigrip bag

**Product characteristics**

Closure type	film
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	full skirt
High or low profile	Low Profile
Number of wells	384
Volume of work	2 µl - 40 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## **Size**

Sample volume	40 µl
Max. volume	40 µl
Diameter	3.2 mm
Width of product	85.5 mm
Height of product	9.65 mm
Length of product	127.8 mm

## **Material & colours**

Product material	Polypropylene (PP)
Colour of product	transparent

## **Purity & certification**

Product category	not regulated
CE certified	no CE mark
Purity standard	PCR Performance Tested
Batched	yes

## **Packaging**

Minimum order qty.	50
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	50
Piece (s) / outer case	50
Piece(s) / pallet	5050
Depth of box	254 mm
Width of box	212 mm
Height of box	165 mm
Depth of case	254 mm
Width of case	212 mm
Height of case	165 mm
Case volume	0.0089 cbm
Weight of product	0.0374 kg
Weight of case	2.18 kg
EAN of inner box	4038917135016
EAN case	4038917067584

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

