



Product description

Order number	72.1984.202
Product description	PCR plate full skirt, 384 well, transparent, Low Profile, 40 µl, PCR Performance Tested, material: PP, 50 piece(s)/minigrip bag

Product characteristics

Closure type	film
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	full skirt
High or low profile	Low Profile
Number of wells	384
Volume of work	2 µl - 40 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Size

Sample volume	40 µl
Max. volume	40 µl
Diameter	3.2 mm
Width of product	85.5 mm
Height of product	9.65 mm
Length of product	127.8 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent

Purity & certification

Product category	no medical device IVD
Purity standard	PCR Performance Tested
Batched	yes

Packaging

Minimum order qty.	50
Type of inner packaging	minigrip bag
Piece(s) / inner box	50
Piece(s) / case	50
Piece(s) / pallet	5050
Case volume	0.0092 cbm
Weight of product	0.0374 kg
Weight of case	2.18 kg
EAN of inner box	4038917135016
EAN case	4038917067584

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

