




**Product Specification**  
**PCR plate full skirt, 384 well, white, Low Profile,**  
**40 µl, PCR Performance Tested, PP**

**Page 1**



**Product description**

Order number	72.1985.202
Product description	PCR plate full skirt, for LightCycler® 480, 384 well, white, Low Profile, 40 µl, PCR Performance Tested, material: PP, optimised for qPCR, 50 piece(s)/minigrip bag

**Product characteristics**

Suitable for	LightCycler® 480
Suitable for qPCR	yes
Base shape	conical
Construction of the plate edges	full skirt
High or low profile	Low Profile
Number of wells	384
Volume of work	2 µl - 40 µl
Centrifugation max (RCF)	3000 x g
Centrifugation temperature	22
Centrifugation time	5 Min
Autoclavable (max. °C)	121 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Size**

Sample volume	40 µl
Max. volume	40 µl
Diameter	3.3 mm
Width of product	85.5 mm
Height of product	9.6 mm
Length of product	127.8 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	white

### **Purity & certification**

Product category	no medical device   IVD
Purity standard	PCR Performance Tested
Batched	yes

### **Packaging**

Minimum order qty.	100
Type of smallest subpackaging	minigrip bag
Piece(s) / inner box	50
Piece (s) / outer case	100
Piece(s) / pallet	7200
Depth of case	279 mm
Width of case	252 mm
Height of case	342 mm
Case volume	0.024 cbm
Weight of product	0.0158 kg
Weight of case	2.15 kg
EAN of inner box	4038917251259
EAN case	4038917251242

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

