



### Product description

Order number 95.1992  
Product description PCR film, for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992), free from DNase/RNase, transparent, suitable for qPCR, self-adhesive, application area: -80 - 120 °C, 100 piece(s)/box

### Product characteristics

Closure type film  
Light protection no  
Suitable for PCR plates (REF 72.1978/.010, 72.1979/.010/.003/.201/.700, 72.1980/.010/.201/.600/.700, 72.1981/.010) & PCR strips (REF 72.985.002/.092/.992)  
Suitable for qPCR yes  
Minimum temperature -80 °C  
Maximum temperature 120 °C

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Colour of product                      transparent

### **Purity & certification**

Product category                      no medical device | IVD  
Purity standard                        free from DNase/RNase  
Batched                                    yes

### **Packaging**

Minimum order qty.                    100  
Type of smallest                        box  
subpackaging  
Piece(s) / inner box                    100  
Piece (s) / outer case                 5000  
Piece(s) / pallet                        90000  
Depth of box                            90 mm  
Width of box                             46 mm  
Height of box                            170 mm  
Depth of case                            380 mm  
Width of case                            380 mm  
Height of case                           380 mm  
Case volume                             0.0549 cbm  
Weight of product                      0.0036 kg  
Weight of case                           20 kg  
EAN of inner box                        4038917940078  
EAN case                                 4038917940085

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

