



### **Product description**

Order number	95.1993
Product description	PCR film, free from DNase/RNase, material: PL, transparent, suitable for qPCR, self-adhesive, application area: -40 - 120 °C, thin material, high optical clarity, 100 piece(s)/bag

### **Product characteristics**

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-40 °C
Maximum temperature	120 °C

### **Size**

Width of product	80 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material                      polyester (PL)  
Colour of product                      transparent

**Purity & certification**

Product category                      not regulated  
CE certified                              no CE mark  
Purity standard                        free from DNase/RNase  
Batched                                  yes

**Packaging**

Minimum order qty.                    100  
Type of smallest                        bag  
subpackaging  
Piece(s) / inner box                    100  
Piece (s) / outer case                 1600  
Piece(s) / pallet                        192000  
Depth of box                            148 mm  
Width of box                            148 mm  
Height of box                            116 mm  
Depth of case                          216 mm  
Width of case                          216 mm  
Height of case                          203 mm  
Case volume                            0.0095 cbm  
Weight of product                      0.002 kg  
Weight of case                         3.74 kg  
EAN of inner box                        4038917303569  
EAN case                                4038917388221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

