



### Product description

Order number	95.1993
Product description	PCR film, free from DNase/RNase, material: PL, transparent, suitable for qPCR, self-adhesive, application area: -40 - 120 °C, thin material, high optical clarity, 100 piece(s)/bag

### Product characteristics

Closure type	film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-40 °C
Maximum temperature	120 °C

### Size

Width of product	80 mm
Length of product	137 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Material & colours**

Product material	polyester (PL)
Colour of product	transparent

**Purity & certification**

Product category	no medical device   IVD
Purity standard	free from DNase/RNase
Batched	yes

**Packaging**

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1600
Piece(s) / pallet	192000
Depth of case	216 mm
Width of case	216 mm
Height of case	203 mm
Case volume	0.0095 cbm
Weight of product	0.002 kg
Weight of case	3.74 kg
EAN of inner box	4038917303569
EAN case	4038917388221

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

