



Product description

Order number	95.1995
Product description	Aluminium film, free from DNase/RNase, material: aluminium, aluminium, self-adhesive, pierceable, application area: -70 - 110 °C, pierceable, light protection for specimens, high chemical resistance, 100 piece(s)/bag

Product characteristics

Closure type	film
Light protection	yes
Suitable for qPCR	no
Minimum temperature	-70 °C
Maximum temperature	110 °C

Size

Width of product	79 mm
Length of product	141 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Material & colours

Product material	aluminium
Colour of product	aluminium

Purity & certification

Product category	no medical device IVD
Purity standard	free from DNase/RNase
Batched	yes

Packaging

Minimum order qty.	100
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1200
Piece(s) / pallet	288000
Case volume	0.0052 cbm
Weight of product	0.0026 kg
Weight of case	3.43 kg
EAN of inner box	4038917088121
EAN case	4038917044295

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

