



Product description

Order number	95.2092
Product description	Film, free of DNase/RNase, transparent, optimised for qPCR, pierceable, application area: -80°C - 120°C

Product characteristics

Type of cap	Film
Light protection	no
Suitable for qPCR	yes
Minimum temperature	-80 °C
Maximum temperature	120 °C

Size

Width of product	82 mm
Length of product	125 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Colour of product transparent

Purity & certification

Product in accordance with MDD or IVDD no medical device | IVD

Purity standard free of DNase/RNase

Lot or serial number yes

Packaging

Minimum order qty. 100

Piece(s) / inner box 100

Piece(s) / case 1200

Piece(s) / pallet 57600

Case volume 0,0142 cbm

Weight of product 0,0031 kg

Weight of case 4,36 kg

EAN of inner box 4038917335676

EAN case 4038917335669

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

