



### **Product description**

Order number	93.1102.001
Product description	Protective cover, material: PC, brown, (LxWxH): 256 x 62 x 72 mm, suitable for stand no 93.1097.100, 1 piece(s)/bag

### **Product characteristics**

Suitable for	stand no 93.1097.100
Autoclavable (max. °C)	121 °C

### **Size**

Width of product	62 mm
Height of product	72 mm
Length of product	256 mm

### **Material & colours**

Product material	Polycarbonate (PC)
Colour of product	brown

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Protective cover, PC, brown, 256 x 62 x 72 mm,**  
**suitable for stand no 93.1097.100**

**Page 2**

**Purity & certification**

Product category	no medical device   IVD
Batched	no

**Packaging**

Minimum order qty.	1
Type of smallest subpackaging	bag
Piece(s) / inner box	1
Piece (s) / outer case	30
Piece(s) / pallet	1080
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.1286 kg
Weight of case	4.61 kg
EAN of inner box	4038917125529
EAN case	4038917020138

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-04-26

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**